

Product Datasheet

Anti-IGF1R/Cy3 (orb983731)

Description	Anti-IGF1R/Cy3
Species/Host	Rabbit
Reactivity	Canine, Human, Mouse
Conjugation	Cy3
Tested Applications	FC, ICC, IF
Immunogen	KLH conjugated synthetic peptide derived from human IGF-1R alpha chain
Target	IGF1R
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	78/150kDa
Uniprot ID	P08069
Dilution Range	Flow-Cyt=1µg/Test,IF=1:100-500,ICC=1:100
Expiration Date	12 months from date of receipt.