
Product Datasheet

**Anti-beta subunit Cholera Toxin/Cy3 [Out of stock]
(orb984363)**

Description

Anti-beta subunit Cholera Toxin/Cy3

Species/Host

Rabbit

Conjugation

Cy3

Immunogen

KLH conjugated synthetic peptide derived from CTXB/beta subunit Cholera Toxin

Target

ctxB

Preservatives

0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.

Form/Appearance

Lyophilized or Liquid

Concentration

1mg/ml

Storage

Store at -20 °C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20°C. When reconstituted in sterile pH 7.4 0.01M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4 °C.

Note

For research use only

Isotype

IgG

Clonality

Polyclonal

MW

12kDa

CAS Number

24730-31-2

Uniprot ID[P01556](#)**Expiration Date**

12 months from date of receipt.