

Product Datasheet

Anti-phospho-GATA2 (Ser401)/Cy3 (orb984672)

Description

Anti-phospho-GATA2 (Ser401)/Cy3

Species/Host

Rabbit

Reactivity

Human

Conjugation

Cy3

Tested

ICC, IF

Applications
Immunogen

KLH conjugated synthesised phosphopeptide derived from human GATA2 around the phosphorylation site of Ser401

Target

GATA2

Preservatives

0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.

Form/Appearance

Lyophilized or Liquid

Concentration

1mg/ml

Storage

Store at -20 °C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20°C. When reconstituted in sterile pH 7.4 0.01M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4 °C.

Note

For research use only

Isotype

IgG

Clonality

Polyclonal

MW

53kDa

Uniprot ID
P23769
Dilution Range

IF=1:100-500,ICC=1:100-500

Expiration Date

12 months from date of receipt.