

Product Datasheet

Myoglobin antibody (For labeling) (orb98724)

Description	Mouse monoclonal to Myoglobin (For labeling)
Species/Host	Mouse
Conjugation	Unconjugated
Target	Myoglobin
Preservatives	PBS pH 7.4
Form/Appearance	Liquid
Storage	-20°C
Note	For research use only
Application notes	Suggested applications: LFA, EIA, CLIA, WB
Isotype	IgG
Clonality	Monoclonal
Source	Ascites
Expiration Date	12 months from date of receipt.