

Product Datasheet

Phospho-SirT1 (Ser27) Antibody (Cy3) (orb988754)

Description	Phospho-SirT1 (Ser27) Antibody (Cy3)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Cy3
Tested Applications	FC
Immunogen	KLH conjugated Synthesised phosphopeptide derived from human SirT1 around the phosphorylation site of Ser27 AS(p-S)PA
Target	SIRT1
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q96EB6
Expiration Date	12 months from date of receipt.