

Product Datasheet

Phospho-Syk (Tyr525 + Tyr526)/Cy3 Antibody (orb988772)

Description	Phospho-Syk (Tyr525 + Tyr526)/Cy3 Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Cy3
Tested Applications	IF
Immunogen	KLH conjugated Synthesised phosphopeptide derived from human Syk around the phosphorylation site of Tyr525/526 EN(p-Y)(p-Y)KA
Target	SYK
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P43405
Expiration Date	12 months from date of receipt.