

Product Datasheet

Anti-Thyroid hormone receptor alpha/Cy3 (orb989817)

Description	Anti-Thyroid hormone receptor alpha/Cy3
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Cy3
Tested Applications	FC, IF
Immunogen	KLH conjugated synthetic peptide derived from human Thyroid hormone receptor alpha
Target	Thyroid hormone receptor alpha
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P10827
Dilution Range	Flow-Cyt=1ug/Test,IF=1:100-500
Expiration Date	12 months from date of receipt.