


Product Datasheet

PLA2G4C antibody (orb99055)

Description	Goat polyclonal to PLA2G4C
Species/Host	Goat
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	PLA2G4C
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Concentration	0.03-0.1 µg/ml
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	This antibody is expected to recognize all reported isoforms (NP_003697.2; NP_001152794.1; NP_001152795.1).WB: Approx 65kDa band observed in Human Brain (Hippocampus), Liver, Pancreas, and Spleen lysates (calculated MW of 62.5kDa according to NP_001152794.1). Recommended concentration: 0.03-0.1 µg/ml
Protein Sequence	NTEVIREYIFDQLR
Clonality	Polyclonal
MW	60.9; 62.5; 59.4
Entrez	8605
Dilution Range	ELISA: 1:128000, WB: 0.03-0.1 µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of Human Pancreas ...