

Product Datasheet

FABP4 antibody (orb99070)

Description	Goat polyclonal antibody to FABP4
Species/Host	Goat
Reactivity	Canine, Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	FABP4
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Concentration	0.03-0.1 µg/ml
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	WB: Approx 14kDa band observed in Human Adipose tissue lysates (calculated MW of 14.7kDa according to NP_001433.1). Recommended concentration: 0.03-0.1 µg/ml
Protein Sequence	SVNGDVITIKSEST
Clonality	Polyclonal
MW	14.7
Entrez	2167
Dilution Range	ELISA: 1:128000, WB: 0.03-0.1 µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of Human Adipose I...