

Product Datasheet

HSP47 Antibody (orb99078)

Description Mouse Anti-Human HSP47 Monoclonal IgG1 Kappa

Species/Host Mouse

Reactivity Human

Conjugation Unconjugated

Tested Applications ICC, WB

Immunogen Human HSP47, full length

Target HSP47

Preservatives PBS pH 7.4, 50% glycerol, 0.09% sodium azide
*Storage buffer may change when conjugated

Concentration 1 mg/ml

Storage -20°C

Note For research use only

Application notes 1 µg/ml of SMC-203 was sufficient for detection of HSP47 in 20 µg of heat shocked HeLa cell lysate by colorimetric immunoblot analysis using Goat anti-mouse IgG:HRP as the secondary antibody.

Isotype IgG1 Kappa

Clonality Monoclonal

MW 47kDa


Uniprot ID [P50454](#)

NCBI [NP_001193943](#)


Entrez [871](#)

Dilution Range WB (1:1000), ICC/IF (1:100)


Expiration Date 12 months from date of receipt.



Western Blot analysis of human epithelia...



Immunofluorescence analysis of heat shoc...



Immunofluorescence analysis of heat shoc...