

Product Datasheet

PDI Antibody (orb99079)

Description

Mouse Anti-Human EndoPDI Monoclonal IgG2a

Species/Host

Mouse

Reactivity

Human

Conjugation

Unconjugated

Tested

ELISA, ICC, IHC, IP, WB

Applications
Immunogen

Synthesized peptide - 12 amino acids long (ADGEDGQDPHSK) corresponding to residues 52 - 63 of the EndoPDI protein

Target

PDI

Preservatives

PBS pH 7.4, 50% glycerol, 0.09% sodium azide
*Storage buffer may change when conjugated

Concentration

1 mg/ml

Storage

-20°C

Note

For research use only

Application notes

1 µg/ml of SMC-204 was sufficient for detection of EndoPDI in 20 µg of heat shocked HeLa cell lysate by colorimetric immunoblot analysis using Goat anti - mouse IgG:HRP as the secondary antibody.

Isotype

IgG2a

Clonality

Monoclonal

MW


48kDa

Uniprot ID
[Q8NBS9](#)
NCBI
[NP_001139021](#)
Entrez
[81567](#)
Dilution Range

WB (1:200)

Expiration Date

12 months from date of receipt.



Immunofluorescence analysis of cervical

...