



Product Datasheet

Dicer Antibody (orb99085)

Description	Mouse Anti-Mouse Dicer Monoclonal IgG1
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, WB
Immunogen	Fusion protein amino acids 1638-1899 of mouse Endoribonuclease Dicer
Target	Dicer
Preservatives	PBS pH 7.4, 50% glycerol, 0.09% sodium azide *Storage buffer may change when conjugated
Concentration	1 mg/ml
Storage	-20°C
Note	For research use only
Application notes	1 µg/ml of SMC-416 was sufficient for detection of Dicer in 20 µg of rat brain lysate by colorimetric immunoblot analysis using Goat anti-mouse IgG:HRP as the secondary antibody.
Isotype	IgG1
Clonality	Monoclonal
MW	215kDa
Uniprot ID	Q8R418
NCBI	NP_683750.2
Entrez	192119
Dilution Range	WB (1:1000), ICC/IF (1:100)
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of neuroblas...



Western Blot analysis of rat brain membr...