

Product Datasheet

Tau protein (phospho-Ser235) antibody (FITC) (orb9913)

Description

Rabbit polyclonal to Tau protein

conjugated to FITC

Species/Host Rabbit

Reactivity Human, Mouse, Rat

Conjugation FITC

Tested Applications IF

Immunogen KLH conjugated Synthesised phosphopeptide derived from human MAPT isoform 2 around the phosphorylation site of Ser235 [PK(P-S)PS]

Target MAPT

Preservatives 0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.

Form/Appearance Lyophilized or Liquid

Concentration - 100 µg (in 100 µl): 1 mg/ml

Storage Store at -20 °C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20°C. When reconstituted in sterile pH 7.4 0.01M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4 °C.

Note For research use only

Isotype IgG

Clonality Polyclonal

MW 83kDa

Uniprot ID **P10636**

Expiration Date 12 months from date of receipt.