


Product Datasheet

Calcium Sensing Receptor (phospho-Thr888) antibody (orb99145)

Description	Calcium Sensing Receptor (phospho-Thr888) antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IF, WB
Note	For research use only
Clonality	Polyclonal
Dilution Range	WB 1:500-1:2000, IF/ICC 1:100-1:500
Expiration Date	12 months from date of receipt.



Western blot analysis of LOVO whole cell...