


Product Datasheet

SOX-9 (phospho-Ser181) antibody (orb99154)

Description	SOX-9 (phospho-Ser181) antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ICC, IF, IHC-P, WB
Note	For research use only
Clonality	Polyclonal
Dilution Range	WB 1:500-1:2000, IHC-P 1:50-1:200, IF/ICC 1:100-1:500
Expiration Date	12 months from date of receipt.



Western blot analysis of NIH-3T3 whole c...