

Product Datasheet

Prothrombin/F2 Antibody (orb99221)

Description	Prothrombin/F2 Antibody
Species/Host	Rabbit
Reactivity	Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	A synthetic peptide corresponding to a sequence at the C-terminus of rat Prothrombin (600-612aa THVFLKRWMQKV).
Form/Appearance	Lyophilized
Concentration	Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.
Storage	Store at -20°C for one year from date of receipt. After reconstitution, at 4°C for one month. It can also be aliquotted and stored frozen at -20°C for six months. Avoid repeated freeze-thaw cycles.
Note	For research use only
Application notes	Tested Species: In-house tested species with positive results. Predicted Species: Species predicted to be fit for the product based on sequence similarities. Other applications have not been tested. Optimal dilutions should be determined by end users. . Add 0.2ml of distilled water will yield a concentration of 500ug/ml.
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	70412 MW
Uniprot ID	P18292
Dilution Range	Western blot, 0.1-0.5µg/ml, Rat
Expiration Date	12 months from date of receipt.



Anti-F2 antibody, orb99221, Western blot...