


Product Datasheet

Mouse IL-12(p40) ELISA Kit (orb99232)

Description	Mouse IL-12(p40) ELISA Kit
Reactivity	Mouse
Range	15.6 pg/ml - 1,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: M23-S335
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, serum and plasma (heparin, EDTA, citrate).
MW	38235 MW
Uniprot ID	P43432
Sensitivity	10 pg/ml
Dilution Range	Dilution ratio of 1:1, concentration in serum and plasma is around 300 pg/ml.
Expiration Date	Please enquire.



Mouse IL-12
(p40)
PicoKine
ELISA Kit
sta...