


Product Datasheet

Human ADAM12 / MCMP ELISA Kit (orb99234)

Description	Human ADAM12 / MCMP ELISA Kit
Reactivity	Human
Range	156 pg/ml - 10,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: CHO; Immunogen sequence: R29-S513
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, serum, plasma (heparin) and urine.
MW	99542 MW
Uniprot ID	O43184
Sensitivity	10 pg/ml
Dilution Range	Dilution ratio of 1:10, concentration in serum and plasma is around 20 ng/ml.
Expiration Date	Please enquire.



Human
ADAM12
PicoKine
ELISA Kit
standard...