


Product Datasheet

Chk1(Phospho-Ser280) antibody (orb99252)

Description	Rabbit polyclonal to Chk1.Serine/threonine-protein kinase which is required for checkpoint-mediated ...
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Note	For research use only
Clonality	Polyclonal
Dilution Range	WB:1:500-1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis on HeLa cell lysat...