


Product Datasheet

EGFR (Phospho-Thr678) antibody (orb99255)

Description	Rabbit polyclonal to EGFR. EGFR is a receptor tyrosine kinase. Receptor for epidermal growth factor (E...
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IF, WB
Note	For research use only
Clonality	Polyclonal
Dilution Range	WB:1:500-1:2000, IF:1:100-1:500
Expiration Date	12 months from date of receipt.



Western blot analysis on 293 treated wit...