


Product Datasheet

beta Amyloid (Phospho-Thr743/668) antibody (orb99260)

Description	Rabbit polyclonal to Amyloid beta A4.APP a cell surface receptor that influences neurite growth, neu...
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IF, IHC-P, WB
Note	For research use only
Clonality	Polyclonal
Dilution Range	WB:1:500-1:2000,IHC-P:1:50-1:200,IF:1:100-1:500
Expiration Date	12 months from date of receipt.



Western blot analysis on HeLa using phos...