


Product Datasheet

Fyn(Phospho-Tyr530) antibody (orb99261)

Description	Rabbit polyclonal to Fyn.Fyn a tyrosine kinase of the Src family.Implicated in the control of cell g...
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Note	For research use only
Clonality	Polyclonal
Dilution Range	WB:1:500-1:2000,IHC-P:1:50-1:200
Expiration Date	12 months from date of receipt.



Western blot analysis on 293 cell lysate...