


Product Datasheet

Tyrosine Hydroxylase (Phospho-Ser31) antibody (orb99262)

Description	Rabbit polyclonal to Tyrosine Hydroxylase. Tyrosine hydroxylase (EC 1.14.16.2) is involved in the con...
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Note	For research use only
Clonality	Polyclonal
Dilution Range	WB:1:500-1:2000,IHC-P:1:50-1:200
Expiration Date	12 months from date of receipt.



Western blot analysis on 293 treated wit...