


Product Datasheet

SP1 (Phospho-Thr453) antibody (orb99263)

Description	Rabbit polyclonal to SP1.SP1 is a transcription factor of the Sp1 C2H2-type zinc-finger protein fami...
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IF, IHC-P, WB
Note	For research use only
Clonality	Polyclonal
Dilution Range	WB:1:500-1:2000,IHC-P:1:50-1:200,IF:1:100-1:500
Expiration Date	12 months from date of receipt.



Western blot analysis on A549 using phos...