


Product Datasheet

Caspase 3 (Phospho-Ser150) antibody (orb99285)

Description	Rabbit polyclonal antibody to Caspase 3 (Phospho-Ser150)
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Note	For research use only
Clonality	Polyclonal
Dilution Range	WB:1:500-1:2000,IHC-P:1:50-1:200
Expiration Date	12 months from date of receipt.



Western blot analysis on Jurkat treated ...