


Product Datasheet

CrkII (Phospho-Tyr221) antibody (orb99286)

Description	Rabbit polyclonal to CrkII. Crk is an adaptor protein with an SH2-SH3-SH3 domain structure. Recruits cyto...
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IF, WB
Note	For research use only
Clonality	Polyclonal
Dilution Range	WB:1:500-1:2000, IF:1:100-1:500
Expiration Date	12 months from date of receipt.



Western blot analysis on COS7 using phos...