


Product Datasheet

GSK3 alpha/ beta (Phospho-Tyr279/216) antibody (orb99288)

Description	Rabbit polyclonal to GSK3 alpha-beta. GSK3A a proline-directed protein kinase of the GSK family. Impl...
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Note	For research use only
Clonality	Polyclonal
Dilution Range	WB:1:500-1:2000,IHC-P:1:50-1:200
Expiration Date	12 months from date of receipt.



Western blot analysis on HeLa using phos...