


Product Datasheet

STAT5B (Phospho-Ser731) antibody (orb99289)

Description	Rabbit polyclonal to STAT5B.The protein encoded by this gene is a member of the STAT family of trans...
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Note	For research use only
Clonality	Polyclonal
Dilution Range	WB:1:500-1:2000,IHC-P:1:50-1:200
Expiration Date	12 months from date of receipt.



Western blot analysis on RAW264.7 treated...