


Product Datasheet

Smad2/3 (Phospho-Thr8) antibody (orb99294)

Description	Rabbit polyclonal to Smad2-3. Smad2 ubiquitously expressed transcription factor phosphorylated and ac...
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Note	For research use only
Clonality	Polyclonal
Dilution Range	WB:1:500-1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis on RAW264.7 using ...