


Product Datasheet

PKC delta(Phospho-Ser645) antibody (orb99298)

Description	Rabbit polyclonal to PKC delta. Protein kinase C (PKC) is a family of serine- and threonine-specific p...
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Note	For research use only
Clonality	Polyclonal
Dilution Range	WB:1:500-1:2000, IHC-P:1:50-1:200
Expiration Date	12 months from date of receipt.



Western blot analysis on NIH-3T3 cell ly...