


Product Datasheet

CaMK2 alpha/ beta/ delta (Phospho-Thr305) (orb99311)

Description	Rabbit polyclonal antibody to CaMK2 alpha/ beta/ delta (Phospho-Thr305)
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Note	For research use only
Clonality	Polyclonal
Dilution Range	WB:1:500-1:2000,IHC-P:1:50-1:200
Expiration Date	12 months from date of receipt.



Western
blot
analysis on
NIH-3T3
using p...