


Product Datasheet

Aurora Kinase antibody (orb99318)

Description	Rabbit polyclonal to Aurora Kinase.The protein encoded by this gene is a cell cycle-regulated kinase...
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Note	For research use only
Clonality	Polyclonal
Dilution Range	WB:1:500-1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis on 293 using Auror...