


Product Datasheet

Synaptotagmin antibody (orb99345)

Description	Rabbit polyclonal to Synaptotagmin.SYT2 May have a regulatory role in the membrane interactions duri...
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Note	For research use only
Clonality	Polyclonal
Dilution Range	WB:1:500-1:2000,IHC-P:1:50-1:200
Expiration Date	12 months from date of receipt.



Western blot analysis on HeLa using Syna...