


Product Datasheet

BIK antibody (orb99386)

Description	Rabbit polyclonal to BIK.The protein encoded by this gene is known to interact with cellular and vir...
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ICC, IF, IHC-P, WB
Note	For research use only
Clonality	Polyclonal
Dilution Range	WB:1:500-1:2000,IHC-P:1:50-1:200,IF:1:100-1:500
Expiration Date	12 months from date of receipt.



Western blot analysis on A549 using BIK ...