

Product Datasheet

LRP5/APC antibody (orb993914)

Description	LRP5/APC antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	APC
Tested Applications	ICC, IF
Immunogen	KLH conjugated synthetic peptide derived from human LRP5
Target	LRP5
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	176kDa
Uniprot ID	075197
Dilution Range	IF=1:100-500,ICC=1:100
Expiration Date	12 months from date of receipt.