



Product Datasheet

GATA4 antibody (orb99436)


Description	GATA4 antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, ICC, IF, WB
Note	For research use only
Clonality	Polyclonal
Dilution Range	WB 1:500-1:2000, IF/ICC 1:100-1:500, ELISA: 1:20000-1:40000
Expiration Date	12 months from date of receipt.



Western blot analysis of COLO205 whole c...



Immunofluorescence analysis of 293 cells...



Immunofluorescence analysis of HepG2 cel...