


Product Datasheet

IkappaB-beta antibody (orb99455)

Description	IkappaB-beta antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, ICC, IF, WB
Note	For research use only
Clonality	Polyclonal
Dilution Range	WB: 1:500-1:2000, IF/ICC: 1:100-1:500, ELISA: 1:20000-1:40000
Expiration Date	12 months from date of receipt.



Western blot analysis of NIH-3T3 cell ly...



Immunofluorescence analysis of HeLa cell...