

AKT (Phospho-Ser473) antibody

Cat#: ORB184255 (CoA)

Description: Rabbit Anti-HSV2 Polyclonal Antibody

Lot No: BR3965

Source: Rabbit

Clonality: Polyclonal

Isotype: IgG

Immunogen: KLH conjugated synthetic peptide derived fom HSV2:11-80/182

Purity: affinity purified by Protein A

Form: Liquid-Buffer: 0.01 TBS(pH 7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol

Intene Use: For Reseach Use o Further manufacture Only. Not for Use in Diagnostic Procedures.

Storage: Store at -20°C. Avoid repeated freeze/thaw cycles.

Biorbyt Ltd.

5 Orwell Furlong, Cowley Road, Cambridge, Cambridgeshire
CB4 0WY, United Kingdom
Email: info@biorbyt.com | Phone: +44 (0)1223 859 353 | Fax: +44(0)1223 280 240

Biorbyt LLC

1100 Corporate Square Drive, Helix Center, Suite 221, St Louis
MO 63132, United States
Email: info@biorbyt.com | Phone: +1 (415)-906-5211 | Fax: +1 (415) 651 8558