

Escherichia coli DH5 alpha antibody

Cat#: orb10594 (CoA)

Name: Escherichia coli DH5 alpha antibody

Description: Rabbit Anti-Escherichia coli DH5 alpha Polyclonal Antibody

Lot No: S4606

Source: Rabbit

Clonality: Polyclonal

Isotype: IgG

Immunogen: whole cell protein of Ecoli DH-5 Alpha

Purity: affinity purified by Protein A

Form: Liquid-Buffer: 0.01M TBS (pH7.4) with 1% BSA, 0.03% Proclin 300 and 50% Glycerol.

Intended use: For Research Use or Further Manufacture Only. Not for Use in Diagnostic Procedures.

Storage: Store at -20°C. Avoid repeated freeze/thaw cycles.

BATCH DATA

IgG/Content Specification: $\geq 1\text{mg/ml}$ (specify method)

IgG/Content Result: 1mg/ml

Antibody Titer Specification: $>10^5$

Antibody Titer Result: $>5.0 \times 10^5$ (ELISA)