

IL18 antibody

Cat#: orb107403 (CoA)

CERTIFICATE OF ANALYSIS

Information

Cat No.: orb107403

Product Name: IL18 Antibody

Lot No.: Q2707

Source: Rabbit Polyconal

Form of Product: Rabbit IgG in pH7.4 PBS, 0.05% NaN₃, 40% Glycerol.

Storage: Store at -20° C/1 year

Quality Analysis

Item Specification Result

ELISA

Dilution: ≥1: 1000 (against immunogen in indirect ELISA)

Dilution: ≥1:2000

Western Blotting (WB)

Dilution: ≥1: 200;

Band Size detected should correspond to the predicted.

Dilution: 1:1350

Predicted band: 22kd

Materials used: 30 µg Hela lysate

Immunohistochemistry (IHC)

Dilution: ≥1: 50;

Cell localization detected should correspond to the predicted.

Dilution: 1:100

Materials used: paraffin-embedded

Human ovarian cancer tissue

Immunofluorescence/ Immunocytochemistry (IF/ICC)

Dilution: ≥1: 50;

Cell localization detected should correspond to the predicted.

Dilution: Not tested

Materials used: /

Caution: This product is for in vitro research use only.