

XIAP antibody

Cat#: orb11576 (CoA)

CERTIFICATE OF ANALYSIS

Description: XIAP Polyclonal Antibody

Lot No. BB10266189

Source: Rabbit

Clonality: Polyclonal

Isotype: IgG

Immunogen: KLH conjugated synthetic peptide derived from mouse XIAP

Purity: affinity purified by Protein A

Form: liquid-Buffer: 0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol

Intended use: For Research Use or Further Manufacture Only. Not for Use in Diagnostic Procedures.

Storage: Shipped at 4C. Store at -20°C. Avoid repeated freeze/thaw cycles.