

Filamin 2 antibody

Cat#: orb156880 (CoA)

Host: Rabbit

Target Protein: FLNC

Lot: M4925

IR: Immunogen Range:251-350/2725

Clonality: Polyclonal

Isotype: IgG

Entrez Gene:[2318](#)

Swiss Prot:[Q14315](#)

Source: KLH conjugated synthetic peptide derived from human Filamin 2:251-350/2725

Purification: affinity purified by Protein A

Storage: 0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol. Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.

Concentration: 1mg/ml

Applications:

WB(1:500-2000)

ELISA(1:5000-10000)

IHC-P(1:100-500)

IHC-F(1:100-500)

ICC(1:100-500)

IF(1:100-500)

Cross Reactive Species:

Human

Mouse

Rat

Dog

Pig

Cow

Horse

Sheep

Chimpanzee.