

Deoxynivalenol antibody

Cat#: orb334823 (CoA)

Catalog #	Source	Reactivity	Applications
orb334823	Mouse	N/A	ELISA
Immunogen	Deoxynivalenol conjugated to KLH		
Purification	Protein G purified		
Specificity	Reacts only with Deoxynivalenol.		
Clonality	Monoclonal		
Form	Liquid in 0.01M Phosphate Buffered Saline, pH 7.2, 50% glycerol, 0.05% Sodium Azide.		
Dilution	ELISA: Use at an assay dependent dilution.		
Gene Symbol			
Storage/Stability	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.		