

Beta-Actin antibody

Cat#: orb499643 (CoA)

CERTIFICATE OF ANALYSIS

Cat. Number: orb499643

Description: beta-Actin (1A2) Monoclonal Antibody

Lot No. BB01056930

Source: Mouse

Clonality: Monoclonal

Isotype: IgG

Immunogen: KLH conjugated synthetic peptide derived from human beta-Actin

Purity: affinity purified by Protein G

Form: 0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol

Intended use: For Research Use or Further Manufacture Only. Not for Use in Diagnostic Procedures.

Storage: Shipped at 4C. Store at -20°C. Avoid repeated freeze/thaw cycles.