

C5b-9 antibody

Cat#: orb499686 (CoA)

Cat No: orb499686

Description: Rabbit Anti-C5b-9 Polyclonal Antibody

Lot No: BA02269454

Source: Rabbit

Clonality: Polyclonal

Isotype: IgG

Immunogen: Purified Human SC5b-9 Complex

Purity: affinity purified by Protein A

From: liquid-Buffer: 0.01M TBS (pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol

Intended use: For Research Use or Further Manufacture Only. Not for Use in Diagnostic Procedures.

Storage: Store at -20°C. Avoid repeated freeze/thaw cycles.

BATCH DATA

IgG/Content Specification: ≥ 1 mg/ml (Specify method)

IgG/Content Result: 1 mg/ml

Antibody Titer Specification: $>10^5$

Antibody Titer Result: $>5.0 \times 10^5$ (ELISA)

Biorbyt Ltd.

5 Orwell Furlong, Cowley Road, Cambridge, Cambridgeshire
CB4 0WY, United Kingdom

Email: info@biorbyt.com | Phone: +44 (0)1223 859 353 | Fax: +44(0)1223 280 240

Biorbyt LLC

1100 Corporate Square Drive, Helix Center, Suite 221, St Louis
MO 63132, United States

Email: info@biorbyt.com | Phone: +1 (415)-906-5211 | Fax: +1 (415) 651 8558