

**CSH2 antibody**

**Cat#: orb52495 (CoA)**

<b>Product Name</b>	CSH2 Antibody
<b>Catalog Number</b>	orb52495
<b>Lot Number</b>	F1214A
<b>Immunogen</b>	Recombinant Human Chorionic somatomammotropin hormone 2 protein (27-217AA)
<b>Host Species</b>	Rabbit
<b>Form</b>	Liquid
<b>Concentration</b>	4 mg/ml
<b>Species Cross-reactivity</b>	Human, Mouse
<b>Purification method</b>	>95%,Protein G purified
<b>Applications</b>	ELISA,WB,IHC;Recommended dilution:WB:1:500-1:5000,IHC:1:20-1:200
<b>Formulation</b>	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
<b>Storage and Stability</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Conclusion</b>	pass