

Quantifoil R 1.2/1.3 on 200 gold mesh

Cat#: orb533752 (User Manual)

Shipping: shipped at ambient temperature

Storage Conditions: store at ambient temperature

Shelf Life: n/a

Applications:

Support film for biological samples in Cryo-EM techniques

Description:

QUANTIFOIL™ R 1.2/1.3 carbon film has circular holes of 1.2 µm diameter with an interspace of 1.3 µm in both dimensions, applied on biocompatible 200 gold mesh. Suitable for magnifications around 50.000x and higher.