

## IL17R antibody

**Cat#: orb5533 (CoA)**

### Certificate of Analysis

**Cat no.:** orb5533

**Description:** Rabbit Anti-IL17RA Polyclonal Antibody

**Lot No.:** BB03025539

**Source:** Rabbit

**Clonality:** Polyclonal

**Isotype:** Rabbit IgG

**Immunogen:** KLH conjugated synthetic peptide derived from human IL-17RA: 201-300/866

**Purification:** Affinity purification

**Form:** Liquid Buffer: 0.01M TBS (pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.

**Intended use:** For Research Use or Further Manufacturer Only. Not for Use in Diagnostic Procedures.

**Storage:** On arrival, Store refrigerated at 2 °C to 8 °C for long term storage product should be frozen in small aliquots between -20 °C to -30 °C. Avoid repeated freeze/thaw cycles.

### Batch Data

**Purity Specification:** ≥90% by SDS-PAGE


**Purity Result:** >95% by SDS-PAGE

**IgG Content Specification:** ≥1.0 mg/ml (specify method)

**IgG Content Result:** 1 mg/ml by BCA Protein Assay Kit (OD562nm)

**Antibody Titer Specification:** >10<sup>5</sup>

**Antibody Titer Result:** >5.0\*10<sup>5</sup> (ELISA)



**Biorbyt Ltd.**

5 Orwell Furlong, Cowley Road, Cambridge, Cambridgeshire  
CB4 0WY, United Kingdom  
Email: info@biorbyt.com | Phone: +44 (0)1223 859 353 | Fax: +44(0)1223 280 240

**Biorbyt LLC**

1100 Corporate Square Drive, Helix Center, Suite 221, St Louis  
MO 63132, United States  
Email: info@biorbyt.com | Phone: +1 (415)-906-5211 | Fax: +1 (415) 651 8558