

## Rabbit Anti BSA antibody

**Cat#: orb572846 (CoA)**

### SPECIFICATION SHEET

**Product Name:** Rabbit Anti BSA

**Catalog #:** orb572846

**Synonyms:** Rabbit anti-BSA(Bovine Serum Albumin)

**Specificity:** Recognizes an epitope on the BSA

**Immunogen:** BSA

**Host Animal:** Rabbit

**Source:** Serum

**Format:** Affinity Purified, Liquid

**Concentration:** 9.3mg/ml **Lot:** LB2209

**Buffer:** 0.01M PBS, pH 7.4

**Preservative:** None

**Description:** Full length protein corresponding to Cow Bovine Serum Albumin aa 1 to the C-terminus. Antisera to bovine albumin was raised by repeated immunisations of rabbits with highly purified bovine albumin. Bovine serum albumin is widely used in biochemical experiments

**Storage & Shelf life:** Store at -20 ° C for 3 years. Avoid multiple freeze/thaw cycles.

#### Biorbyt Ltd.

5 Orwell Furlong, Cowley Road, Cambridge, Cambridgeshire  
CB4 0WY, United Kingdom

Email: [info@biorbyt.com](mailto:info@biorbyt.com) | Phone: +44 (0)1223 859 353 | Fax: +44(0)1223 280 240

#### Biorbyt LLC

Suite 103, 369 Pine Street, San Francisco  
California 94104, United States

Email: [info@biorbyt.com](mailto:info@biorbyt.com) | Phone: +1 (415)-906-5211 | Fax: +1 (415) 651 8558