

Diphtheria Toxoid antibody (HRP)

Cat#: orb593632 (CoA)

Catalog Number orb593632

Product Name Mouse Anti-Diphtheria Toxoid MAb, HRP-conjugated

Batch number KB1368

Clonity Monoclonal

Size 100ug, 500ug

Immunogen Diphtheria Toxoid

Species Mouse

Tag HRP

Buffer 0.1M PBS, PH 8.0, with protein stabilizer

Purification Dialysis after Labeling

Background The filamentous haemagglutinin adhesin (FHA) is a large, filamentous protein that serves as a dominant attachment factor for adherence to host ciliated epithelial cells of the respiratory tract, called respiratory epithelium. It is associated with biofilm formation and possesses at least four binding domains which can bind to different cell receptors on the epithelial cell surface.

Storage Store at 4°C for 6 months. It is recommended to aliquote to ~10ul per vial and store at -20°C for long term storage. Avoid repeated freezing and thawing cycles.

Biorbyt Ltd.

5 Orwell Furlong, Cowley Road, Cambridge, Cambridgeshire
CB4 0WY, United Kingdom

Email: info@biorbyt.com | Phone: +44 (0)1223 859 353 | Fax: +44(0)1223 280 240

Biorbyt LLC

Suite 103, 369 Pine Street, San Francisco
California 94104, United States

Email: info@biorbyt.com | Phone: +1 (415)-906-5211 | Fax: +1 (415) 651 8558